



IHE Integration Statement			
Vendor	Product Name	Version	Date
ScImage, Inc.	PicomEnterprise	3.3.x	December 14, 2016
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Scheduled Workflow (SWF)	Image Manager/Image Archive		
	Performed Procedure Step Manager		
	Image Display		
	Evidence Creator		
Patient Information Reconciliation (PIR)	Image Manager/Image Archive		
	Performed Procedure Step Manager		
Consistent Presentation of Images (CPI)	Image Manager/Image Archive		
	Image Display		
	Evidence Creator		
Access to Radiology Information (ARI)	Report Reader	Multiple Sources	
	Report Repository		
	External Report Repository Access		
	Image Display	Multiple Sources	
	Image Manager/Image Archive		
Key Image Note (KIN)	Evidence Creator		
	Image Manager/Image Archive		
	Image Display		
Simple Image and Numeric Report (SINR)	Report Creator		
	Report Manager		
	Report Repository		
	Report Reader		
Teaching File and Clinical Trial Export (TCE)	Export Selector		
	Export Manager		
Cross-Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source	PDF Report	
	Imaging Document Consumer	Retrieve Images Retrieve Evidence Documents	
Mammography Image (MAMMO)	Image Manager/Image Archive		
	Image Display		
	Evidence Creator		
Image Object Change Management (IOCM)	Image Manager/Image Archive		
	Image Display		
Cardiac Catheterization Workflow (CATH)	Image Manager/Image Archive	Cardiac CATH Intermittently Connected Modality	
	Performed Procedure Step Manager		
	Image Display		
	Evidence Creator		
Echocardiography Workflow (ECHO)	Image Manager/Image Archive	Echocardiography Intermittently Connected Modality	
	Performed Procedure Step Manager		
	Image Display	Stress Echo	
	Evidence Creator		
Consistent Time (CT)	Time Client	SNTP	
Internet address for vendor's IHE information: www.scimage.com/support/ScImage_IHE_Statement.pdf			



Links to Standards Conformance Statements for the Implementation		
HL7	http://www.scimage.com/support-documentation/HL7.pdf	
DICOM	www.scimage.com/support/PicomEnterprise3.3DCS.pdf	
Links to general information on IHE		
In North America: www.ihe.net	In Europe: www.ihe-europe.net	In Japan: www.ihe-j.org